

Alpha Liner

Alpha Original MP Liner

The Alpha Original MP Liner is the powerhouse core of the Alpha Liner Family. Since its introduction in 1996 the Alpha Original Liner has brought comfort and protection to amputees worldwide.

Made from thermoplastic elastomer, the Alpha gel gently adheres to the skin to protect against abrasion and breakdown while also proving to be extremely skin-friendly. The Alpha gel is covered in a layer of durable fabric that helps to extend the life of the liner and makes it very easy for amputees to slide into their sockets. The Alpha Original MP Liner is heat moldable, allowing prosthetists to easily customize the liner within their facilities.



Product No.	Description
ALC-5063-E	Alpha Original MP Cushion Liner, Small
ALC-5060-E	Alpha Original MP Cushion Liner, Medium
ALC-5064-E	Alpha Original MP Cushion Liner, Medium Plus
ALC-5066-E	Alpha Original MP Cushion Liner, Large
ALC-5068-E	Alpha Original MP Cushion Liner, Large Plus
ALC-5067-E	Alpha Original MP Cushion Liner, Extra Large
ALL-5063-E	Alpha Original MP Locking Liner, Small
ALL-5060-E	Alpha Original MP Locking Liner, Medium
ALL-5064-E	Alpha Original MP Locking Liner, Medium Plus
ALL-5066-E	Alpha Original MP Locking Liner, Large
ALL-5068-E	Alpha Original MP Locking Liner, Large Plus
ALL-5067-E	Alpha Original MP Locking Liner, Extra Large

**Locking or Cushion
Liner - Same Price!**

SPECIFICATIONS

MODEL	Original
TYPE	Locking or Cushion
STYLES	Uniform
COLORS	Original (Green and Gray)
THICKNESS	6 mm
SIZES	Small, Medium, Medium Plus, Large, Large Plus, Extra Large
U.S. CODE SUGGESTION	L-5679 Cushion Liner (bill once per liner) L-5673 Locking Liner (bill once per liner)

For locking liners, Ohio Willow Wood strongly recommends applying Loc-Tite Removable Threadlocker 242 (or equivalent) to the pin threads prior to installing the pin. Failure to ensure that the locking pin is secure in the Locking Liner could result in the pin coming out of the liner, resulting in the potential for injury to the amputee.